

Product Datasheet

ERV3 antibody (orb125089)

Description Goat Polyclonal to ERV3. NP_001007254.2, PMID: 15994821 and GeneID: 2086.

Species/Host Goat

Reactivity Human

Conjugation Unconjugated

Tested Applications ELISA, IHC, WB

Target ERV3

Preservatives Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.

Storage Aliquot and store at -20°C. Minimize freezing and thawing.

Note For research use only

Application notes WB: This product has been successfully used in WB on Human: Nahas MR et al, Br J Haematol. 2019 May;185(4):679-690. PMID: 30828801.. IHC: Paraffin embedded Human Placenta. Recommended concentration: 3.75µg/ml. This product has been successfully used in IHC on Human: Frese S et al, PLoS One. 2015 Jul 8;10(7). PMID: 26154387.

Protein Sequence PQDNFTLTASS


Clonality Polyclonal

MW 67.9

Entrez **2086**

Dilution Range ELISA: 1:4000, IHC-P: 3.75 µg/ml

Expiration Date 12 months from date of receipt.



Immunohistochemical staining of Human Placenta