


Product Datasheet

IL5 Antibody (Biotin) (orb1250905)

Description	IL5 Antibody (Biotin)
Species/Host	Rat
Reactivity	Human
Conjugation	Biotin
Tested Applications	FC
Immunogen	COS-expressed human IL-5
Target	IL5
Preservatives	Supplied in PBS/NaN ₃
Form/Appearance	Liquid
Concentration	0.5 mg/mL
Storage	Store vial at 2-8°C
Note	For research use only
Application notes	IL-5 Antibody [JES1-39D10] for use in flow cytometry, ELISA, and immunohistochemistry / immunocytochemistry assays.
Isotype	IgG2a
Clonality	Monoclonal
Uniprot ID	P05113
NCBI	P05113
Dilution Range	IL-5 Antibody [JES1-39D10] for use in flow cytometry, ELISA, and immunohistochemistry / immunocytochemistry assays.
Expiration Date	12 months from date of receipt.



Standard curve generated with Rat Anti-H...