


Product Datasheet

IFNG Antibody (orb1250919)

Description	IFNG Antibody
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC
Immunogen	E. coli-expressed human IFN-gamma
Target	IFNG
Preservatives	Supplied in borate buffered saline, pH 8.2; No preservatives or amine-containing buffer salts added
Form/Appearance	Liquid
Concentration	0.5 mg/mL
Storage	Store vial at 2-8°C
Note	For research use only
Application notes	IFN-gamma Antibody [A35] for use in ELISA assays.
Isotype	IgG1
Clonality	Monoclonal
Uniprot ID	P01579
NCBI	P01579
Dilution Range	IFN-gamma Antibody [A35] for use in ELISA assays.
Expiration Date	12 months from date of receipt.



Standard curve generated with Mouse Anti...