


Product Datasheet

IL10 Antibody (Biotin) (orb1250924)

Description	IL10 Antibody (Biotin)
Species/Host	Rat
Reactivity	Human, Primate, Virus
Conjugation	Biotin
Tested Applications	Cell Separation, FC, IHC-Fr, IP
Immunogen	COS-expressed human IL-10
Target	IL10
Preservatives	Supplied in PBS/NaN ₃
Form/Appearance	Liquid
Concentration	0.5 mg/mL
Storage	Store vial at 2-8°C
Note	For research use only
Application notes	IL-10 Antibody [JES3-12G8] for use in immunohistochemistry, ELISA, ELISpot, western blotting, and multiplex assays.
Isotype	IgG2a
Clonality	Monoclonal
Uniprot ID	P22301
NCBI	P22301
Dilution Range	IL-10 Antibody [JES3-12G8] for use in immunohistochemistry, ELISA, ELISpot, western blotting, and multiplex assays.
Expiration Date	12 months from date of receipt.



Standard curve generated with Rat Anti-H...