

Product Datasheet

FAS Antibody (orb1250929)

Description	FAS Antibody
Species/Host	Mouse
Reactivity	Human, Rat
Conjugation	Unconjugated
Tested Applications	FC
Immunogen	Human CD95 transfected L cells
Target	FAS
Preservatives	Supplied in PBS
Form/Appearance	Liquid
Concentration	0.5 mg/mL
Storage	Store vial at -20°C
Note	For research use only
Application notes	CD95 Antibody [DX3] for use in flow cytometry, ELISA, immunohistochemistry assays.
Isotype	IgG2a
Clonality	Monoclonal
Uniprot ID	P25445
NCBI	P25445
Dilution Range	CD95 Antibody [DX3] for use in flow cytometry, ELISA, immunohistochemistry assays.
Expiration Date	12 months from date of receipt.