

Product Datasheet

FAS Antibody (orb1250930)

Description	FAS Antibody
Species/Host	Mouse
Reactivity	Human, Monkey, Primate
Conjugation	Unconjugated
Tested Applications	FC, IP
Immunogen	Human CD95 transfected L cells
Target	FAS
Preservatives	Supplied in PBS
Form/Appearance	Liquid
Concentration	0.5 mg/mL
Storage	Store vial at -20°C
Note	For research use only
Application notes	CD95 Antibody [DX2] for use in flow cytometry, immunohistochemistry / immunocytochemistry, immunoprecipitation, ELISA, and apoptosis assays.
Isotype	IgG1
Clonality	Monoclonal
Uniprot ID	P25445
NCBI	P25445
Dilution Range	CD95 Antibody [DX2] for use in flow cytometry, immunohistochemistry / immunocytochemistry, immunoprecipitation, ELISA, and apoptosis assays.
Expiration Date	12 months from date of receipt.