


Product Datasheet

CD1C Antibody (FITC) (orb1250943)

Description	CD1C Antibody (FITC)
Species/Host	Mouse
Reactivity	Gallus
Conjugation	FITC
Tested Applications	FC
Immunogen	Lymphocytes from the bursa of Fabricius of outbred chickens
Target	CD1C
Preservatives	Supplied in PBS/NaN ₃
Form/Appearance	Liquid
Concentration	0.5 mg/mL
Storage	Store vial at 2-8°C
Note	For research use only
Application notes	CD1.1 Antibody [CB3] for use in flow cytometry and immunohistochemistry / immunocytochemistry assays.
Isotype	IgG1
Clonality	Monoclonal
Dilution Range	CD1.1 Antibody [CB3] for use in flow cytometry and immunohistochemistry / immunocytochemistry assays.
Expiration Date	12 months from date of receipt.



Chicken
peripheral
blood
mononuclear
cel...