


Product Datasheet

EIF3E antibody (orb125095)

Description	Goat polyclonal antibody to EIF3E
Species/Host	Goat
Reactivity	Bovine, Canine, Human, Porcine
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Target	EIF3E
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	WB: Approx 50kDa band observed in lysates of cell line Daudi (calculated MW of 52.2kDa according to NP_001559.1). Recommended concentration: 0.03-0.1µg/ml.
Protein Sequence	KLNQNSRSEAPN
Clonality	Polyclonal
MW	52.2
Entrez	3646
Dilution Range	ELISA: 1:2000, WB: 0.03-0.1 µg/ml
Expiration Date	12 months from date of receipt.



Western blot analysis of Daudi lysate us...