


Product Datasheet

CD44 Antibody (Biotin) (orb1250954)

Description	CD44 Antibody (Biotin)
Species/Host	Mouse
Reactivity	Aves, Gallus
Conjugation	Biotin
Tested Applications	Blocking, COSTIM, FC, ICC, IHC, IP
Target	CD44
Preservatives	Supplied in PBS/NaN ₃
Form/Appearance	Liquid
Concentration	0.5 mg/mL
Storage	Store vial at 2-8°C
Note	For research use only
Application notes	CD44 Antibody [AV6] for use in flow cytometry and immunohistochemistry assays.
Isotype	IgG1
Clonality	Monoclonal
Dilution Range	CD44 Antibody [AV6] for use in flow cytometry and immunohistochemistry assays.
Expiration Date	12 months from date of receipt.



Chicken
peripheral
blood
mononuclear
cel...