


Product Datasheet

CNTF antibody (orb125096)

Description	Goat polyclonal antibody to CNTF
Species/Host	Goat
Reactivity	Canine, Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Target	CNTF
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	WB: Approx 17+28kDa bands observed in Human, Mouse and Rat Brain lysates (calculated MW of 22.9kDa according to Human NP_000605.1 and 22.6kDa according to Mouse NP_740756.1). This molecular weight is routinely observed by other sources. Recommended concentrat
Protein Sequence	SIHDLRFISSHQ
Clonality	Polyclonal
MW	22.9
Entrez	1270
Dilution Range	ELISA: 1:2000, WB: 0.5-2 µg/ml
Expiration Date	12 months from date of receipt.



Western blot analysis of Mouse Brain lys...