


Product Datasheet

IDH1 antibody (orb125098)

Description	Goat polyclonal antibody to IDH1
Species/Host	Goat
Reactivity	Canine, Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Target	IDH1 (aa153-164)
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	WB: Approx 48kDa band observed in lysates of cell line HepG2 and of Mouse Liver (calculated MW of 46.7kDa according to NP_005887.2). Recommended concentration: 0.3-1µg/ml.
Protein Sequence	EITYTPSDGTQK
Clonality	Polyclonal
MW	46.7
Entrez	3417
Dilution Range	ELISA: 1:1000, WB: 0.3-1 µg/ml
Expiration Date	12 months from date of receipt.



Western blot analysis of HepG2 lysate us...