


Product Datasheet

IgA Antibody [A-1] (Biotin) (orb1250987)

Description	IgA Antibody [A-1] (Biotin)
Species/Host	Mouse
Reactivity	Gallus
Conjugation	Biotin
Tested Applications	FC, ICC, IHC-Fr, IP
Preservatives	Supplied in PBS/NaN3
Form/Appearance	Liquid
Concentration	0.5 mg/mL
Storage	Store vial at 2-8°C
Note	For research use only
Application notes	IgA Antibody [A-1] for use in flow cytometry, immunohistochemistry / immunocytochemistry, ELISA, and immunoprecipitation assays.
Isotype	IgG2b
Clonality	Monoclonal
Dilution Range	IgA Antibody [A-1] for use in flow cytometry, immunohistochemistry / immunocytochemistry, ELISA, and immunoprecipitation assays.
Expiration Date	12 months from date of receipt.



Chicken intestine tissue was stained with...