


Product Datasheet

KLHDC8B antibody (orb125100)

Description	Goat Polyclonal to KLHDC8B. NP_775817.1, PMID: 20107318 and GeneID: 200942.
Species/Host	Goat
Reactivity	Bovine, Canine, Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Target	KLHDC8B
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	WB: Approx 38kDa band observed in lysates of cell lines Daudi, Jurkat and NIH3T3 (calculated MW of 37.7kDa according to NP_775817.1). Recommended concentration: 0.5-2µg/ml. An additional band of unknown identity was also consistently observed at 22kDa. This b
Protein Sequence	TKLPRSLRMRDK
Clonality	Polyclonal
MW	37.7
Entrez	200942
Dilution Range	ELISA: 1:64000, WB: 0.5-2 µg/ml
Expiration Date	12 months from date of receipt.



Western blot analysis of Daudi lysate us...