


Product Datasheet

TCRab/Vb1 Antibody [TCR-2] (FITC) (orb1251017)

Description	TCRab/Vb1 Antibody [TCR-2] (FITC)
Species/Host	Mouse
Reactivity	Aves, Gallus
Conjugation	FITC
Tested Applications	CMCD, FC, ICC, IP
Immunogen	White Leghorn chick thymocytes and blood mononuclear cells
Preservatives	Supplied in PBS/NaN3
Form/Appearance	Liquid
Concentration	0.5 mg/mL
Storage	Store vial at 2-8°C
Note	For research use only
Application notes	TCRab/Vb1 Antibody [TCR-2] for use in flow cytometry and immunohistochemistry assays.
Isotype	IgG1
Clonality	Monoclonal
Dilution Range	TCRab/Vb1 Antibody [TCR-2] for use in flow cytometry and immunohistochemistry assays.
Expiration Date	12 months from date of receipt.



Chicken
peripheral
blood
mononuclear
cel...