



Product Datasheet

CD3D Antibody (orb1251030)


Description	CD3D Antibody
Species/Host	Mouse
Reactivity	Aves, Gallus
Conjugation	Unconjugated
Tested Applications	FC, IHC-Fr
Immunogen	Chicken thymocytes and Ig- blood mononuclear cells
Target	CD3D
Preservatives	Supplied in borate buffered saline, pH 8.2; No preservatives or amine-containing buffer salts added
Form/Appearance	Liquid
Concentration	0.5 mg/mL
Storage	Store vial at 2-8°C
Note	For research use only
Application notes	CD3 Antibody [CT-3] for use in flow cytometry, immunohistochemistry, immunoprecipitation, and stimulation assays.
Isotype	IgG1
Clonality	Monoclonal
Dilution Range	CD3 Antibody [CT-3] for use in flow cytometry, immunohistochemistry, immunoprecipitation, and stimulation assays.
Expiration Date	12 months from date of receipt.



Paraffin
embedded
chicken
spleen tissue
...



Frozen vIBDV
infected
bursal tissue
sect...



Frozen
chicken ovary
section was
stained...