

Product Datasheet

CD1 Antibody (FITC) (orb1251069)

Description	CD1 Antibody (FITC)
Species/Host	Mouse
Reactivity	Porcine
Conjugation	FITC
Tested Applications	FLISA
Immunogen	Fresh dd miniature swine thymocytes
Target	CD1
Preservatives	Supplied in PBS/NaN ₃
Form/Appearance	Liquid
Concentration	0.5 mg/mL
Storage	Store vial at 2-8°C
Note	For research use only
Application notes	CD1 Antibody [76-7-4] for use in flow cytometry and immunohistochemistry / immunocytochemistry assays.
Isotype	IgG2a
Clonality	Monoclonal
Uniprot ID	Q9XS72
NCBI	Q9XS72
Dilution Range	CD1 Antibody [76-7-4] for use in flow cytometry and immunohistochemistry / immunocytochemistry assays.
Expiration Date	12 months from date of receipt.