


Product Datasheet

Cd86 Antibody (FITC) (orb1251093)

Description	Cd86 Antibody (FITC)
Species/Host	Rat
Reactivity	Mouse
Conjugation	FITC
Tested Applications	FC, ICC, IHC-Fr, IHC-P, IP
Immunogen	Mouse B cell lymphoma cell line 5C2
Target	Cd86
Preservatives	Supplied in PBS/NaN3
Form/Appearance	Liquid
Concentration	0.5 mg/mL
Storage	Store vial at 2-8°C
Note	For research use only
Application notes	CD86 Antibody [2D10] for use in flow cytometry and immunohistochemistry assays.
Isotype	IgG2b
Clonality	Monoclonal
Uniprot ID	P42082
NCBI	P42082
Dilution Range	CD86 Antibody [2D10] for use in flow cytometry and immunohistochemistry assays.
Expiration Date	12 months from date of receipt.



DBA/2 mouse mastocytoma cell line P815
a...