



Product Datasheet

SOX2 antibody (orb12511)

Description	Rabbit polyclonal antibody to SOX2
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IF, IHC, WB
Immunogen	Recombinant fusion protein of human SOX2(NP_003097.1).
Preservatives	PBS with 0.05% proclin300,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	35KDa
Uniprot ID	P48431
Dilution Range	WB: 1:500-1:2000
Expiration Date	12 months from date of receipt.




Immunofluorescence analysis of HeLa cell...



Immunofluorescence analysis of C6 cells

...



Western blot analysis of extracts of var...