

Product Datasheet

H-2Kd/H-2Dd Antibody [34-1-2S] (FITC) (orb1251106)

Description	H-2Kd/H-2Dd Antibody [34-1-2S] (FITC)
Species/Host	Mouse
Reactivity	Mouse
Conjugation	FITC
Tested Applications	FC
Immunogen	BDF1 mouse splenocytes
Target	H2-D1; H2-K1
Preservatives	Supplied in PBS/NaN ₃
Form/Appearance	Liquid
Concentration	0.5 mg/mL
Storage	Store vial at 2-8°C
Note	For research use only
Application notes	H-2Kd/H-2Dd Antibody [34-1-2S] for use in flow cytometry and immunocytochemistry assays.
Isotype	IgG2a
Clonality	Monoclonal
Dilution Range	H-2Kd/H-2Dd Antibody [34-1-2S] for use in flow cytometry and immunocytochemistry assays.
Expiration Date	12 months from date of receipt.