

Product Datasheet

H-2Kd/H-2Dd Antibody [34-1-2S] (orb1251107)

Description	H-2Kd/H-2Dd Antibody [34-1-2S]
Species/Host	Mouse
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	FC, IP
Immunogen	BDF1 mouse splenocytes
Target	H2-D1; H2-K1
Preservatives	Supplied in borate buffered saline, pH 8.2; No preservatives or amine-containing buffer salts added
Form/Appearance	Liquid
Concentration	0.5 mg/mL
Storage	Store vial at 2-8°C
Note	For research use only
Application notes	H-2Kd/H-2Dd Antibody [34-1-2S] for use in flow cytometry, immunocytochemistry, electron microscopy, immunoprecipitation, CMCD, blocking, and adhesion assays.
Isotype	IgG2a
Clonality	Monoclonal
Dilution Range	H-2Kd/H-2Dd Antibody [34-1-2S] for use in flow cytometry, immunocytochemistry, electron microscopy, immunoprecipitation, CMCD, blocking, and adhesion assays.
Expiration Date	12 months from date of receipt.