


Product Datasheet

NOX1 antibody (orb125114)


Description	Goat polyclonal antibody to NOX1
Species/Host	Goat
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Target	NOX1
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	WB: Approx 70-72kDa band observed in lysates of cell lines A549, CaCo-2, HeLa and U937, (calculated MW of 64.9kDa according to NP_008983.2;). This molecular weight is routinely observed by other sources and was successfully blocked by incubation with the immu
Protein Sequence	DKATDIVTGLKQK
Clonality	Polyclonal
MW	64.9; 59.0
Entrez	27035
Dilution Range	ELISA: 1:128000, WB: 0.1-0.3 µg/ml, IHC-P: 7.5 µg/ml
Expiration Date	12 months from date of receipt.



Western blot analysis of colorectal adenocarcinoma tissue.



Immunohistochemical staining of Human Kidney tissue.



Immunohistochemical staining of Human Kidney tissue.