


Product Datasheet

POLDIP2 antibody (orb125116)

Description	Goat Polyclonal to POLDIP2. NP_056399.1, PMID: 20158880 and GeneID: 26073.
Species/Host	Goat
Reactivity	Bovine, Canine, Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Target	POLDIP2
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	WB: Approx 37kDa band observed in Rat Vascular Smooth Muscle Cells lysates (calculated MW of 15.2kDa according to Rat NP_001099286.1). The band diminished upon siRNA knock down. Recommended concentration: 1-3µg/ml.
Protein Sequence	RFERPDGSHFDVR
Clonality	Polyclonal
MW	42
Entrez	26073
Dilution Range	ELISA: 1:32000, WB: 1-3 µg/ml
Expiration Date	12 months from date of receipt.



Western blot analysis of Rat VSMC lysate...