



Product Datasheet

POLDIP2 antibody (orb125117)


Description	Goat polyclonal antibody to POLDIP2
Species/Host	Goat
Reactivity	Bovine, Canine, Human, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Target	POLDIP2
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	WB: Approx 37kDa band observed in lysates of cell lines K562 and MCF7 and of Mouse Liver (calculated MW of 42.0kDa according to NP_056399.1). Approx 37kDa band observed in Rat Vascular Smooth Muscle Cells lysates. The band diminished upon siRNA knock down. Reco
Protein Sequence	KTHTYYQVLIDARDC
Clonality	Polyclonal
MW	42
Entrez	26073
Dilution Range	ELISA: 1:8000, WB: 1-3 µg/ml
Expiration Date	12 months from date of receipt.



Western blot analysis of staining of K562...



Western blot analysis of staining of Rat...



Western blot analysis of staining of Mou...