


Product Datasheet

Arl6 antibody (orb125119)

Description	Goat Polyclonal to Arl6. NP_062639.3, PMID: 20207729.
Species/Host	Goat
Reactivity	Bovine, Canine, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC
Target	Arl6 / BBS3 (mouse)
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	WB: Preliminary experiments gave an approx 33kDa band in Mouse Brain lysates after 1µg/ml antibody staining. Please note that currently we cannot find an explanation in the literature for the band we observe given the calculated size of 21.0kDa according to N. IHC: Paraffin embedded Human Lung. Recommended concentration: 5µg/ml.
Protein Sequence	KPSNAQSQDIVPT
Clonality	Polyclonal
MW	21
Dilution Range	ELISA: 1:16000, WB: 1 µg/ml, IHC-P: 5 µg/ml
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of Human Lu...