

Product Datasheet

Ctla4 Antibody (Biotin) (orb1251208)

Description	Ctla4 Antibody (Biotin)
Species/Host	Hamster
Reactivity	Mouse
Conjugation	Biotin
Tested Applications	FC
Immunogen	Extracellular portion of murine CTLA-4 fused to a murine IgG2a
Target	Ctla4
Preservatives	Supplied in PBS/NaN ₃
Form/Appearance	Liquid
Concentration	0.5 mg/mL
Storage	Store vial at 2-8°C
Note	For research use only
Application notes	CD152 Antibody [1B8] for use in flow cytometry and ELISA assays.
Isotype	IgG
Clonality	Monoclonal
Uniprot ID	P09793
NCBI	P09793
Dilution Range	CD152 Antibody [1B8] for use in flow cytometry and ELISA assays.
Expiration Date	12 months from date of receipt.