


Product Datasheet

ODF2 antibody (orb125122)

Description	Goat Polyclonal to ODF2. NP_002531.3; NP_702915.1, PMID: 19386263 and GeneID: 4957.
Species/Host	Goat
Reactivity	Canine, Human, Mouse, Porcine, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Target	ODF2
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	WB: Approx 85kDa band observed in C57/B6 mouse CTL lysates (calculated MW of 95.4kDa according to Mouse NP_001106684.1). Recommended concentration: 1-3µg/ml.
Protein Sequence	REKHQASQKENKQ
Clonality	Polyclonal
MW	93.3; 73.3
Entrez	4957
Dilution Range	ELISA: 1:128000, WB: 1-3 µg/ml
Expiration Date	12 months from date of receipt.



Western blot analysis of mouse CTL lysat...