


Product Datasheet

Thy1 Antibody (FITC) (orb1251233)

Description	Thy1 Antibody (FITC)
Species/Host	Rat
Reactivity	Mouse
Conjugation	FITC
Tested Applications	FC, IHC-Fr, WB
Immunogen	Mouse thymus or spleen
Target	Thy1
Preservatives	Supplied in PBS/NaN3
Form/Appearance	Liquid
Concentration	0.5 mg/mL
Storage	Store vial at 2-8°C
Note	For research use only
Application notes	CD90.2 Antibody [30-H12] for use in flow cytometry and immunohistochemistry / immunocytochemistry assays.
Isotype	IgG2b
Clonality	Monoclonal
Uniprot ID	P01831
NCBI	P01831
Dilution Range	CD90.2 Antibody [30-H12] for use in flow cytometry and immunohistochemistry / immunocytochemistry assays.
Expiration Date	12 months from date of receipt.



BALB/c
mouse
splenocytes
were stained
wi...