

Product Datasheet

Muc19 antibody (orb125125)

Description	Goat polyclonal antibody to Muc19
Species/Host	Goat
Reactivity	Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA
Target	Mucin 19 / Apomucin
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	WB: Not yet tested - our routinely used western blotting protocol does not allow detection of proteins as large as the calculated size of 694kDa according to NP_997126.2. Therefore we cannot recommend an optimal concentration and the antibody is an aspiring p
Protein Sequence	DTFRDCHSKVDPTTY
Clonality	Polyclonal
MW	694
Dilution Range	ELISA: 1:128000
Expiration Date	12 months from date of receipt.