

Product Datasheet

Bpifa2e antibody (orb125127)

Description	Goat Polyclonal to Bpifa2e. NP_032979.1, PMID: 19499239.
Species/Host	Goat
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	ELISA
Target	Parotid secretory protein (mouse)
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	WB: Not yet tested. At this stage we are dependent on researchers in the field for further characterization of this product. Therefore we cannot recommend an optimal concentration and the antibody is an aspiring product. We would appreciate any feedback from
Protein Sequence	DLINSLSIKTNAQ
Clonality	Polyclonal
MW	24.8
Dilution Range	ELISA: 1:8000
Expiration Date	12 months from date of receipt.