



Product Datasheet

Psp antibody (orb125128)

Description	Goat Polyclonal to Psp. NP_032979.1, PMID: 19499239.
Species/Host	Goat
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, ICC, WB
Target	parotid secretory protein (mouse aa146-160)
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	WB: Approx 25kDa band observed in Mouse mandibular glands from E17 onwards (calculated MW of 24.8kDa according to NP_032979.1). Recommended concentration: 1-3µg/ml. Data provided by Everest Grant winner Melinda Larsen, State University of New York, Albany, NY. Notes: Immunocytochemistry: A subset of murine submandibular gland cells (E18) show specific staining in secretory granules. Recommended concentration 4ug/ml. Data provided by Everest Grant winner Melinda Larsen, State University of New York, Al
Protein Sequence	KTNAQTGLPEVTIGK
Clonality	Polyclonal
MW	24.8
Dilution Range	ELISA: 1:128000, WB: 1-3 µg/ml, IHC-P: 4ug/ml
Expiration Date	12 months from date of receipt.



Western blot analysis of Mouse Mandibular glands.



Immunofluorescence analysis of mouse Submandibular glands.