


Product Datasheet

Cd40Ig Antibody (FITC) (orb1251285)

Description	Cd40Ig Antibody (FITC)
Species/Host	Hamster
Reactivity	Mouse
Conjugation	FITC
Tested Applications	FC
Immunogen	Activated mouse Th1 clone D1.6
Target	Cd40Ig
Preservatives	Supplied in PBS/NaN ₃
Form/Appearance	Liquid
Concentration	0.5 mg/mL
Storage	Store vial at 2-8°C
Note	For research use only
Application notes	CD154 Antibody [MR1] for use in flow cytometry and immunohistochemistry assays.
Isotype	IgG
Clonality	Monoclonal
Uniprot ID	P27548
NCBI	P27548
Dilution Range	CD154 Antibody [MR1] for use in flow cytometry and immunohistochemistry assays.
Expiration Date	12 months from date of receipt.



PMA and Ionomycin-stimulated BALB/c mous...