


Product Datasheet

Dennd4c antibody (orb125129)

Description	Goat Polyclonal to Dennd4c. NP_908976.1; NP_001074483.1, PMID: 14610273.
Species/Host	Goat
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Target	Dennd4c
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	WB: Approx 210kDa band observed in Mouse Spleen lysates (calculated MW of 211kDa according to NP_908976.1). Recommended concentration: 0.5-2µg/ml.
Protein Sequence	RFKPQSPYRAYKDRS
Clonality	Polyclonal
MW	211; 179
Dilution Range	ELISA: 1:8000, WB: 0.5-2 µg/ml
Expiration Date	12 months from date of receipt.



Western blot analysis of Mouse Spleen ly...