


## Product Datasheet

### CD16/32 Antibody [93] (Biotin) (orb1251300)

|                            |                                                                                                        |
|----------------------------|--------------------------------------------------------------------------------------------------------|
| <b>Description</b>         | CD16/32 Antibody [93] (Biotin)                                                                         |
| <b>Species/Host</b>        | Rat                                                                                                    |
| <b>Reactivity</b>          | Mouse                                                                                                  |
| <b>Conjugation</b>         | Biotin                                                                                                 |
| <b>Tested Applications</b> | FC                                                                                                     |
| <b>Immunogen</b>           | VH81X Tg B cells and T220 stromal cell line                                                            |
| <b>Target</b>              | Fcgr2b; Fcgr3                                                                                          |
| <b>Preservatives</b>       | Supplied in PBS/NaN3                                                                                   |
| <b>Form/Appearance</b>     | Liquid                                                                                                 |
| <b>Concentration</b>       | 0.5 mg/mL                                                                                              |
| <b>Storage</b>             | Store vial at 2-8°C                                                                                    |
| <b>Note</b>                | For research use only                                                                                  |
| <b>Application notes</b>   | CD16/32 Antibody [93] for use in flow cytometry, immunohistochemistry, and immunoprecipitation assays. |
| <b>Isotype</b>             | IgG2a                                                                                                  |
| <b>Clonality</b>           | Monoclonal                                                                                             |
| <b>Uniprot ID</b>          | <b>P08101, P08508</b>                                                                                  |
| <b>NCBI</b>                | <b>P08101, P08508</b>                                                                                  |
| <b>Dilution Range</b>      | CD16/32 Antibody [93] for use in flow cytometry, immunohistochemistry, and immunoprecipitation assays. |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                                                        |



BALB/c  
mouse  
splenocytes  
were stained  
wi...