


Product Datasheet

Cd28 Antibody (FITC) (orb1251313)

Description	Cd28 Antibody (FITC)
Species/Host	Hamster
Reactivity	Mouse
Conjugation	FITC
Tested Applications	FC, ICC, IHC-Fr, IHC-P
Immunogen	C57BL/6N mouse T cell lymphoma EL4 cell line
Target	Cd28
Preservatives	Supplied in PBS/NaN ₃
Form/Appearance	Liquid
Concentration	0.5 mg/mL
Storage	Store vial at 2-8°C
Note	For research use only
Application notes	CD28 Antibody [PV-1] for use in flow cytometry assays.
Isotype	IgG
Clonality	Monoclonal
Uniprot ID	P31041
NCBI	P31041
Dilution Range	CD28 Antibody [PV-1] for use in flow cytometry assays.
Expiration Date	12 months from date of receipt.



BALB/c mouse thymocytes were stained with...