


Product Datasheet

Ascl3 antibody (orb125133)

Description	Goat polyclonal antibody to Ascl3
Species/Host	Goat
Reactivity	Bovine, Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Target	Ascl3 (mouse)
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	WB: Approx 23kDa band observed in Mouse Liver lysates (calculated MW of 20.2kDa according to NP_064435.1). Recommended concentration: 0.1-0.3µg/ml. Primary incubation was 1 hour.. IHC: Paraffin embedded Human Brain (Cortex). Recommended concentration: 3.75µg/ml.
Protein Sequence	HLPEDYLEKRLSK
Clonality	Polyclonal
MW	20.2
Entrez	56676
Dilution Range	ELISA: 1:8000, WB: 0.1-0.3 µg/ml, IHC-P: 3.75 µg/ml
Expiration Date	12 months from date of receipt.



Western blot analysis of Mouse Liver lys...