

Product Datasheet

Itgb2 Antibody (FITC) (orb1251348)

Description	Itgb2 Antibody (FITC)
Species/Host	Rat
Reactivity	Mouse
Conjugation	FITC
Tested Applications	FC, ICC, IHC-Fr
Immunogen	Cell membrane glycoproteins derived from mouse BW5147 T cell lymphoma cells
Target	Itgb2
Preservatives	Supplied in PBS/NaN ₃
Form/Appearance	Liquid
Concentration	0.5 mg/mL
Storage	Store vial at 2-8°C
Note	For research use only
Application notes	CD18 Antibody [C71/16] for use in flow cytometry, immunohistochemistry / immunocytochemistry, and western blotting assays.
Isotype	IgG2a
Clonality	Monoclonal
Uniprot ID	P11835
NCBI	P11835
Dilution Range	CD18 Antibody [C71/16] for use in flow cytometry, immunohistochemistry / immunocytochemistry, and western blotting assays.
Expiration Date	12 months from date of receipt.