



Product Datasheet

Itgam Antibody (FITC) (orb1251356)


Description	Itgam Antibody (FITC)
Species/Host	Rat
Reactivity	Mouse
Conjugation	FITC
Tested Applications	ELISA
Immunogen	Peritoneal macrophages from B6D2 hybrid mice
Target	Itgam
Preservatives	Supplied in PBS/NaN3
Form/Appearance	Liquid
Concentration	0.5 mg/mL
Storage	Store vial at 2-8°C
Note	For research use only
Application notes	CD11b Antibody [3A33] for use in flow cytometry and immunohistochemistry / immunocytochemistry assays.
Isotype	IgG2a
Clonality	Monoclonal
Uniprot ID	P05555
NCBI	P05555
Dilution Range	CD11b Antibody [3A33] for use in flow cytometry and immunohistochemistry / immunocytochemistry assays.
Expiration Date	12 months from date of receipt.



BALB/c mouse splenocytes were stained wi...



Frozen thermally burned mouse tissue sec...



Frozen thermally burned mouse tissue sec...