


Product Datasheet

Vcam1 Antibody (FITC) (orb1251401)

Description	Vcam1 Antibody (FITC)
Species/Host	Rat
Reactivity	Mouse
Conjugation	FITC
Tested Applications	ELISA, WB
Immunogen	BALB/3T3 and +/+2.4 cells
Target	Vcam1
Preservatives	Supplied in PBS/NaN3
Form/Appearance	Liquid
Concentration	0.5 mg/mL
Storage	Store vial at 2-8°C
Note	For research use only
Application notes	CD106 Antibody [M/K-2] for use in flow cytometry, immunohistochemistry, and western blotting assays.
Isotype	IgG1
Clonality	Monoclonal
Uniprot ID	P29533
NCBI	P29533
Dilution Range	CD106 Antibody [M/K-2] for use in flow cytometry, immunohistochemistry, and western blotting assays.
Expiration Date	12 months from date of receipt.



C57BL/6
mouse bone
marrow cells
were sta...