


Product Datasheet

Sap30L antibody (orb125141)

Description	Goat Polyclonal to Sap30L. NP_001074637.1, PMID: 19015240 and GenElD: 79685.
Species/Host	Goat
Reactivity	Canine, Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Target	Sap30l (mouse)
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	WB: Approx 18kDa band observed in Rat Spleen lysates (calculated MW of 20.8kDa according to Rat XP_340804.4). Recommended concentration: 0.3-1µg/ml.
Protein Sequence	RLDQKSEGSQLE
Clonality	Polyclonal
MW	20.7
Entrez	79685
Dilution Range	ELISA: 1:32000, WB: 0.3-1 µg/ml
Expiration Date	12 months from date of receipt.



Western blot analysis of Rat Pancreas ly...