


Product Datasheet


Type V Collagen Antibody (orb1251416)

Description	Type V Collagen Antibody
Species/Host	Goat
Reactivity	Canine, Equine, Fish, Gallus, Guinea pig, Human, Mammal, Mouse, Porcine, Rabbit, Rat
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	Type V collagen
Preservatives	Supplied in borate buffered saline, pH 8.2; No preservatives or amine-containing buffer salts added
Form/Appearance	Liquid
Concentration	0.4 mg/mL
Storage	Store vial at 2-8°C
Note	For research use only
Application notes	Type V Collagen Antibody for use in ELISA, immunohistochemistry / immunocytochemistry, electron microscopy, and western blotting assays.
Isotype	IgG
Clonality	Polyclonal
Dilution Range	Type V Collagen Antibody for use in ELISA, immunohistochemistry / immunocytochemistry, electron microscopy, and western blotting assays.
Expiration Date	12 months from date of receipt.




Human gastric cancer tissue was stained

...



Frozen D-arginine treated Mdx mouse quad...



Paraffin embedded mouse inguinal mammary...