



Product Datasheet

Type III Collagen Antibody (orb1251420)


Description	Type III Collagen Antibody
Species/Host	Goat
Reactivity	Bovine, Canine, Equine, Gallus, Guinea pig, Hamster, Human, Mammal, Mouse, Porcine, Primate, Rabbit, Rat, Rodent
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	Type III collagen
Preservatives	Supplied in borate buffered saline, pH 8.2; No preservatives or amine-containing buffer salts added
Form/Appearance	Liquid
Concentration	0.4 mg/mL
Storage	Store vial at 2-8°C
Note	For research use only
Application notes	Type III Collagen Antibody for use in ELISA, immunohistochemistry / immunocytochemistry, electron microscopy, flow cytometry, and western blotting assays.
Isotype	IgG
Clonality	Polyclonal
Dilution Range	Type III Collagen Antibody for use in ELISA, immunohistochemistry / immunocytochemistry, electron microscopy, flow cytometry, and western blotting assays.
Expiration Date	12 months from date of receipt.



Human gastric cancer tissue was stained ...



Frozen osteopontin-null mouse muscle sec...



Rat kidney section was stained with Goat...