


Product Datasheet

Cyp2d5 antibody (orb125143)

Description	Goat Polyclonal to Cyp2d5. NP_775426.1, PMID: 2819073.
Species/Host	Goat
Reactivity	Rat
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Target	Cyp2d5 / Cyp2d10 (rat)
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	WB: Approx 65kDa band observed in Rat Spleen lysates (calculated MW of 57.1kDa according to NP_775426.1). Recommended concentration: 1-3µg/ml.
Protein Sequence	CDAFTAQNGRSINPK
Clonality	Polyclonal
MW	57.1
Dilution Range	ELISA: 1:32000, WB: 1-3 µg/ml
Expiration Date	12 months from date of receipt.



Western blot analysis of Rat Spleen lysa...