

Product Datasheet

SGCD antibody (orb125144)

Description

Goat Polyclonal to SGCD. NP_000328.2;

Species/Host

GeneID: 6444.

Goat

Reactivity

Bovine, Canine, Human, Mouse, Porcine, Rat

Conjugation

Unconjugated

Tested

ELISA, IHC, WB

Applications
Target

Delta-Sarcoglycan (aa245-257)

Preservatives

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.

Storage

Aliquot and store at -20°C. Minimize freezing and thawing.

Note

For research use only

Application notes

WB: Approx 37kDa band observed in Human and Rat Heart lysates (calculated MW of 32.2kDa according to NP_000328.2). The observed molecular weight corresponds to the glycosylated form. Recommended concentration: 0.1-0.3µg/ml.. IHC: In paraffin embedded Human Skeletal Muscle shows sarcolemmar staining of the fibres. Recommended concentration, 3-5µg/ml.

Protein Sequence

RLPHGSYPTGTR

Clonality

Polyclonal

MW

32.2; 28.8; 32.1

Entrez


6444

Dilution Range


ELISA: 1:16000, WB: 0.1-0.3 µg/ml, IHC-P: 3-5µg/ml

Expiration Date

12 months from date of receipt.

Western blot analysis of human heart lys...



Immunohistochemical analysis of formalin...