


Product Datasheet

BID Antibody (orb1251448)

Description	BID Antibody
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC
Immunogen	Recombinant Bid
Target	BID
Preservatives	Supplied in borate buffered saline, pH 8.2; No preservatives or amine-containing buffer salts added
Form/Appearance	Liquid
Concentration	0.5 mg/mL
Storage	Store vial at 2-8°C
Note	For research use only
Application notes	Bid Antibody [SB50a] for use in ELISA, western blotting, and immunoprecipitation assays.
Isotype	IgG2b
Clonality	Monoclonal
Uniprot ID	P55957
NCBI	P55957
Dilution Range	Bid Antibody [SB50a] for use in ELISA, western blotting, and immunoprecipitation assays.
Expiration Date	12 months from date of receipt.



Total cell lysates from Jurkat cells wer...