


Product Datasheet

IL10 Antibody (FITC) (orb1251460)

Description	IL10 Antibody (FITC)
Species/Host	Rat
Reactivity	Mouse
Conjugation	FITC
Tested Applications	FC
Immunogen	E. coli-expressed mouse IL-10
Target	IL10
Preservatives	Supplied in PBS/NaN ₃
Form/Appearance	Liquid
Concentration	0.1 mg/mL
Storage	Store vial at 2-8°C
Note	For research use only
Application notes	IL-10 Antibody [JES5-16E3] for use in flow cytometry and immunohistochemistry / immunocytochemistry assays.
Isotype	IgG2b
Clonality	Monoclonal
Uniprot ID	P18893
NCBI	P18893
Dilution Range	IL-10 Antibody [JES5-16E3] for use in flow cytometry and immunohistochemistry / immunocytochemistry assays.
Expiration Date	12 months from date of receipt.



PMA and
Ionomycin-
stimulated
BALB/c
mous...