


Product Datasheet

STK38 antibody (orb125147)

Description	Goat polyclonal antibody to STK38
Species/Host	Goat
Reactivity	Bovine, Canine, Human, Porcine, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Target	NDR1 / STK38
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	WB: Approx 55kDa band observed in Human and Rat Spleen lysates (calculated MW of 54.2kDa according to NP_009202.1). Recommended concentration: 0.5-2µg/ml.
Protein Sequence	KPTVATSNHPET
Clonality	Polyclonal
MW	54.2
Entrez	11329
Dilution Range	ELISA: 1:32000, WB: 0.5-2 µg/ml
Expiration Date	12 months from date of receipt.



Western blot analysis of Human Spleen ly...