


Product Datasheet

BCL2L1 Antibody (PE) (orb1251472)

Description	BCL2L1 Antibody (PE)
Species/Host	Mouse
Reactivity	Human, Mouse, Primate, Rat
Conjugation	PE
Tested Applications	FC, ICC, IHC-Fr, IHC-P
Immunogen	Recombinant human Bcl-xS
Target	BCL2L1
Preservatives	Supplied in PBS/NaN3 and a stabilizing agent
Form/Appearance	Liquid
Concentration	0.1 mg/mL
Storage	Store vial at 2-8°C
Note	For research use only
Application notes	Bcl-xL Antibody [7B2.5] for use in flow cytometry assays.
Isotype	IgG3
Clonality	Monoclonal
Uniprot ID	Q07817
NCBI	Q07817
Dilution Range	Bcl-xL Antibody [7B2.5] for use in flow cytometry assays.
Expiration Date	12 months from date of receipt.



T cell
leukemia
cell line
Jurkat was
int...