

Product Datasheet

FAS Antibody (orb1251484)

Description

FAS Antibody

Species/Host

Mouse

Reactivity

Human, Monkey, Primate

Conjugation

Unconjugated

Tested

FC, IP

Applications
Immunogen

Human CD95 transfected L cells

Target

FAS

Preservatives

Supplied in borate buffered saline, pH 8.2; No preservatives or amine-containing buffer salts added

Form/Appearance

Liquid

Concentration

0.1 mg/mL

Storage

Store vial at 2-8°C

Note

For research use only

Application notes

CD95 Antibody [DX2] for use in flow cytometry, immunohistochemistry / immunocytochemistry, immunoprecipitation, ELISA, and apoptosis assays.

Isotype

IgG1

Clonality

Monoclonal

Uniprot ID
P25445
NCBI
P25445
Dilution Range

CD95 Antibody [DX2] for use in flow cytometry, immunohistochemistry / immunocytochemistry, immunoprecipitation, ELISA, and apoptosis assays.

Expiration Date

12 months from date of receipt.