

Product Datasheet

ITGAL Antibody (orb1251526)

Description	ITGAL Antibody
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC
Immunogen	Fibronectin purified monocytes
Target	ITGAL
Preservatives	Supplied in borate buffered saline, pH 8.2; No preservatives or amine-containing buffer salts added
Form/Appearance	Liquid
Concentration	0.1 mg/mL
Storage	Store vial at 2-8°C
Note	For research use only
Application notes	CD11a Antibody [38] for use in flow cytometry, immunohistochemistry, immunoprecipitation, ELISA, blocking, and adhesion assays.
Isotype	IgG2a
Clonality	Monoclonal
Uniprot ID	P20701
NCBI	P20701
Dilution Range	CD11a Antibody [38] for use in flow cytometry, immunohistochemistry, immunoprecipitation, ELISA, blocking, and adhesion assays.
Expiration Date	12 months from date of receipt.