

Product Datasheet

IgY Antibody [G-1] (PE) (orb1251546)

Description	IgY Antibody [G-1] (PE)
Species/Host	Mouse
Reactivity	Aves, Gallus
Conjugation	PE
Tested Applications	CMCD, FC, IHC-Fr
Preservatives	Supplied in PBS/NaN ₃ and a stabilizing agent
Form/Appearance	Liquid
Concentration	0.1 mg/mL
Storage	Store vial at 2-8°C
Note	For research use only
Application notes	IgY Antibody [G-1] for use in flow cytometry assays.
Isotype	IgG1
Clonality	Monoclonal
Dilution Range	IgY Antibody [G-1] for use in flow cytometry assays.
Expiration Date	12 months from date of receipt.

